# SOPHIA GENETICS

Developed and tested by genomic experts, the SOPHiA GENETICS<sup>™</sup> Pharmacogenomics Community Panel combines a capture-based target enrichment kit with the analytical capabilities and advanced features of the SOPHiA DDM<sup>™</sup> Platform.

## **Main Features**

The SOPHiA GENETICS<sup>®</sup> Pharmacogenomics Community Panel covers coding and non-coding regions of 41 genes associated with pharmacogenomics. Probe design is highly optimized to guarantee a high on-target reads percentage and coverage uniformity even in GC-rich regions, including the first exon. A specialized analytical module covers CNVs, star alleles, and promoters/UTRs in *CYP2D6*, and the pharmacogenomics-specific output files report the genotype status of all genes in the panel.

Gene Panel	Variants Called	Recommendations	Wet Lab
ABCB1, ABCC2, ABCC4, ABCG2, ACE, ADRB1,	SNVs	Starting material	Day 1:
ADRB2, BCHE, CDA, COMT, CYP1A2, CYP2B6,	Indels	200 ng DNA	Library Preparation
CYP2C19, CYP2C8, CYP2C9, CYP2D6, CYP3A4, CYP3A5, CYP4F2, DPYD, DRD2, G6PD, HTR2A, ITPA, MTHFR, NAT2, NR1I2, NR1I3, NUDT15,	CNVs Star alleles ( <i>CYP2D6</i> only)	Sample type Blood	Day 2: Capture and Sequencing
OPRM1, P2RY12, POR, RYR1, SLC22A1, SLC01B1, TPMT, TYMS, UGT1A1, UGT1A4, UGT2B7, VKORC1	. ,	Samples per run for > 50x coverage depth 48 for Illumina MiSeq® V3 32 for Illumina MiSeq® V2 48 for Illumina MiniSeq™ High Output Kit 16 for Illumina MiniSeq™ Mid Output Kit	Hands-on library preparation time: 8 hours

## **Analytical Performance**

The SOPHiA DDM<sup>™</sup> Platform analyzes complex NGS data by detecting and annotating multiple types of genomic variants in all the genes of the panel.

#### Analysis time from FASTQ: 2 hours

	Observed	Lower 95% Cl
Sensitivity	100%	96.34%
Accuracy	100%	99.99%
Precision	100%	96.34%
Specificity	100%	99.99%
Repeatability	99.98%	99.97%
Reproducibility	99.98%	99.97%
Average on-target rate	67%	
Coverage uniformity	99.78%	
Average % of target region with depth >200x	99.83%	

## One Simple Intuitive Platform: Beyond Analytics

## Accelerated assessment and reporting of genomic variants

Dedicated features in SOPHiA DDM<sup>™</sup> reduce the complexity of determining the significance of genomic variants and facilitate the interpretation process, thus reducing turnaround time:

- Variant annotation Using International guidelines such as ACMG and information from curated databases such as ClinVar
- Variant filtering Virtual Panels and Cascading Filters for speedy screening of relevant variants
- Variant status report Pharmacogenomics-specific output files with genotype status

## Global support at every step

We offer local support anywhere in the world. Our dedicated bioinformaticians help save time and resources, ensuring fast resolution of workflow disruptions. In addition, our Set Up Program provides assistance with assay set up for a fast and worry-free transition to routine testing.

## Secure and unlimited data storage

Access to the SOPHiA DDM" Platform is restricted to registered users only. The Platform provides unlimited and unrestricted storage, while keeping data safe by applying the highest industrial standards of encryption in compliance with your local data security policies.

## Access to the SOPHiA GENETICS community

In the SOPHiA DDM" Platform, experts from hundreds of healthcare institutions interpret their results and flag the pathogenicity level of variants according to their knowledge and experience. This highly valuable information feeds the variant knowledge base and is anonymously and safely shared among the members of the community.

Product code: CS2532ILLRGLY10

1. Excluding CYP2D6 analytical module. Analytical performance metrics were calculated from 2 runs with 12 samples on an Illumina MiSeq<sup>®</sup> instrument, and 3 runs with 12 samples on an Illumina MiniSeq<sup>™</sup> instrument.

For Research Use Only. Not for Use in Diagnostic Procedures.

SOPHIA GENETICS does not intend to use or display any company trademarks or trade names to imply a relationship with, or endorsement or sponsorship of SOPHIA GENETICS or its products by any company.